

Introducing...

The **NEW CondenSyn**

Distillation Adapter

New application tool developed in partnership with the University of East Anglia

The NEW CondenSyn Distillation Adapter allows chemists to use the CondenSyn waterless air condenser in a vertical orientation for distillation purposes.

This new application represents an extension to the capabilities and potential water savings delivered by the CondenSyn in laboratories worldwide.

Our customers are concerned about both the environmental impact and indeed running costs of running a research facility where a condenser is considered an essential tool. The water-free alternative, CondenSyn, not only reduces water waste dramatically but is easy to clean, fast to set up, simple to swap between sizes, and has a handy non-roll feature to stop accidents when left on a bench.

Paired with the new Distillation Adapter , CondenSyn provides unbeatable value for the planet and for your budget!



Tests were completed using an IKA RCT basic with PT1000 temperature probe, DrySyn Scholar heat block, 250 ml round bottom flask, cone thread adapter fitted with thermometer and 100 ml round bottom collection flask. For experiments using the conventional apparatus a still head, water condenser, and vented bend were used. For the CondenSyn experiments the distillation attachment and a CondenSyn replaced these items, pictured right.







SUSTAINABLE LABS

Effective. Sustainable. Safe. Brilliant.

sales@quarkglass.com Tel: 800-955-0376



Introducing...

The **NEW CondenSyn**

Distillation Adapter

Test Data:

Experiment	Duration (mins)	Condenser Type	% Loss
Clove oil distillation	90	CondenSyn	0.88
		Water Condenser	1.40
Ethanol removal (50:50 aqueous base/ethanol)	60	CondenSyn	0.63
		Water Condenser	1.17
Methanol removal (90:10 ethan-1,2-diol/methanol	60	CondenSyn	0.62
		Water Condenser	0.72
Cyclopentadiene distillation	60	CondenSyn	0.64
		Water Condenser	0.51
Pyrrole vacuum distillation	30	CondenSyn	2.01
		Water Condenser	1.66

synt

Choose your CondenSyn:

US-M-200-1420 200mm with 14/20 Joint US-M-200-1922 200mm with 19/22 Joint US-M-200-2440 200mm with 24/40 Joint

Then add the Distillation Adapter

US-C-DIST-24

CondenSyn Distillation Adapter with 24/40 fitttings



SUSTAINABLE LABS

Effective. Sustainable. Safe. Brilliant.

sales@quarkglass.com Tel: 800-955-0376