

## The DrySyn Parallel Synthesis Kit & Starter Kit

**No more oil baths or clamps needed!**

The DrySyn Parallel Synthesis Kit allows the scientist to interchangeably heat and stir multiple reactions in any standard size tube or vial using any standard hotplate.

**Complete kit includes:** MULTI-E base plate to fit any hotplate stirrer, heat resistant handle and 1 pack of three standard reaction vial inserts in one size.

**Starter kit includes:** MULTI-E Base with heat resistant handle, and your choice of any three standard reaction vial inserts in any size from those detailed below.

**Remember to specify which size inserts you require when ordering!**



### Popular purchase options

**Complete kit part number: ADS3-PS-(size)**

**Starter kit part number: ADS3-PSS-(sizes)**

Hole size (Ø mm)	Holes per insert	Suited for:
11.80 Ø	9	2 mL HPLC sample vials
15.00 Ø	5	1 dram vials
16.20 Ø	4	CEM microwave tubes
17.40 Ø	4	Biotage 2-5 mL microwave tubes & 2 dram vials
20.20 Ø	4	Robot tubes
20.40 Ø	4	24 mm boiling tubes
25.75 Ø	4	Wheaton reactor vials
27.65 Ø	4	20 mL Scintillation vials
28.20 Ø	4	Biotage 10-20 mL microwave tubes

### Also Available

**Part Number: ADS3-R-(size)**

Pack of 3 x standard reaction vial inserts in **one size** from those detailed here.

**Part Number: ADS3-R-A-(sizes)**

Pack of 3 x standard reaction vial inserts in **any size** from those detailed here.

**Part Number: ADS3-RT-17,20**

Pack of 3 x tapered **17.20 mm** inserts for Biotage 0.5-2 mL 17.00 mm microwave vials

(Not currently available as part of the Parallel Synthesis kits)



**Also available custom manufactured to your required size!**

To enquire about all your available options contact Quark Glass  
800-955-0376, sales@quarkglass.com